EN – ENGLISH – Instructions for use

- Do not use if package is damaged*
- Terminology and package contents**
- A. Fixing handle
 - B. Needle pusher
 - C. Distal tag pusher
 - D. Proximal tag pusher
 - E. Safety pin
 - F. Brake
 - G. Needle and tags (in the dispenser)
 - H. Dispenser attachment

- Technical information**
- Needle diameter – 1.75 mm / 5,25 Fr
 - Needle rigid length – 17 mm
 - Total needle length (excl. handle) – 1 200 mm
 - Tag diameter – 1.1 mm / 3,3 Fr
 - Tag length – 10 mm
 - Suture type – EP1.5 (USP 4-0) braided polyester
 - Maximum duration of use – 4 hours

Intended use

TAPES is a suturing device using a flexible endoscope. For use in natural orifices.

Contraindications

Any contraindication to endoscopy; gastric ulcer; perforations; severe gastritis; bleeding; oesophagitis; colitis; hiatal hernia; women who are pregnant or trying to become pregnant; addiction (drug or alcohol); bleeding disorders, rectitis; varices; portal hypertension in the stomach; fibrotic sites

Warnings

In order to ensure the safe use of TAPES, it is essential to read this paragraph carefully.

- Before using TAPES, make sure that the package is properly sealed and includes all the components shown in the picture and that they are in perfect condition. If any component is missing or damaged, do not use TAPES.
- TAPES is a device for gastroenterology. It can only be used in natural orifices.
- TAPES can only be used by qualified medical staff (a gastroenterologist or paramedical staff under surveillance by a gastroenterologist) who has received the necessary training and certification to use it.
- Always keep at least one additional TAPES in stock before starting a procedure.
- This document describes the use of TAPES, not the clinical procedures that can be performed with it.
- TAPES is single-use and cannot be reused, resterilised or reloaded.
- TAPES can only be used on patients under general anaesthesia.
- TAPES can only be used with adult triangulation platforms. See "Triangulation platform compatibility" for more details.
- If TAPES is placed outside a clean zone, it cannot be used.
- If one of the components of TAPES becomes detached and falls into a body cavity of the patient, retrieve it with endoscopic forceps if it is judged dangerous for the patient's health.
- If the triangulation platform does not allow large enough plications to be formed, the suture may not hold.

Precautions

- Never make sudden movements, especially during insertion and removal. This could hurt the patient.
- The needle must not be apparent during the introduction of the triangulation platform. This could hurt the patient.
- TAPES must always be manipulated under visual surveillance when the needle is outside the triangulation platform, in the cavity to be treated.
- Always wear appropriate protective clothing, as blood and other tissues carry an infection risk. Suitable protective clothing includes protective goggles, a mask and waterproof clothes and gloves.
- Do not use force to insert TAPES into the triangulation system. If there is resistance, remove it and investigate.
- When multiple tag chains are used, avoid sutures crossing each other.
- Do not insert TAPES into a bent channel. Always introduce TAPES into a channel that is as straight as possible.
- TAPES is a 5.25 Fr device. Use it only in a suitable channel (up to 7 Fr). To use it in a larger channel, use a catheter to avoid guiding issues.
- Always keep the handle firmly in the hand to avoid undesired movements.

Adverse reactions

The following adverse reactions have been observed during endoscopic interventions involving needles and may happen after the use of TAPES:

- Nausea, vomiting, cramps
- Bleeding
- Minor peritonitis
- Minor bacteraemia

Procedure

The procedure described here is a check-list for an operation by an experienced user.

1. Check that the TAPES package is not damaged and is perfectly sealed. If any anomaly is detected, do not use TAPES; send it back to Endo Tools Therapeutics.
2. Insert a compatible triangulation platform at the site to be treated.
3. Unpack TAPES and check visually that it is complete and in good condition. Check that the tags are in place and the handle has not been actuated. Remove TAPES from the dispenser by turning the fixing handle.
4. Check that the needle pusher slides in the fixing handle and check the brake.
5. Gently introduce TAPES into the triangulation platform, using lubricating gel if necessary. Do not use force. Check that the angling of the platform allows insertion.

- EN – Instructions for use
- FR – Mode d’emploi
- NL – Gebruikshandleiding
- DE – Bedienungsanleitung
- IT – Istruzioni per l’uso

6. If necessary, fix TAPES to the triangulation platform with the connector.

7. Actuate the triangulation platform to position the TAPES needle correctly.

8. Check visually that the TAPES needle slides correctly out of the platform.

9. Pierce the first tissue by pushing the handle.

10. Remove the safety pin and release the distal tag using the distal tag pusher.

11. Pull back the distal tag pusher flush with the proximal tag pusher.

12. Pull the needle back into the platform by pulling the handle.

13. Move the platform back to loosen the suture if necessary.

14. Pierce the second tissue.

15. Release the proximal tag using the proximal tag pusher. Check that the entire suture has been released. Make back and forth movements with the tag pusher if necessary.

16. Pull back the proximal tag pusher to its rest position.

17. Using the handle, pull back the needle into the platform.

18. Straighten the triangulation platform to allow needle removal. Always keep visual control.

19. Take TAPES out of the platform. Do not use force.

20. With rat tooth forceps in a catheter, tighten the tag chain by pulling on the suture loop. Tighten until the tags exert suitable pressure on the tissues.

21. With a cutting tool, cut the free end of the suture. Make sure you do not cut the suture between tags.

End of procedure

1. Check that all the components are present (excluding tags and suture).
2. Dispose of TAPES in an appropriate trash bin (for sharp hospital waste).

Emergency procedure

If TAPES gets stuck in the triangulation platform or tissue, use the following procedure:

If TAPES is stuck in tissue:

1. Release the proximal tag to release the needle from the tissues.
2. If necessary, cut the tag or suture with an endoscopic guillotine or scissors to release the needle.
3. If necessary, retrieve the tags and suture with endoscopic forceps.

If TAPES is stuck in the triangulation platform:

1. Completely remove the distal and proximal tag handles.
2. Check that the angle of the platform allows TAPES to slide.
3. If TAPES is still stuck, the needle does not slide out of the operating channel and the triangulation platform allows it, remove everything together (platform and TAPES) and investigate the reasons for the blockage.
4. If the needle is out of the channel or the platform does not allow for combined removal, cut the needle from its tube e.g. with argon plasma.
5. Retrieve the needle with the endoscope and foreign body forceps.
6. Remove TAPES from the operating channel. If TAPES is damaged during the operation and a part stays in the stomach, take it out of the patient.

Storage

TAPES must be stored in its sterile package, in a clean and dry place, away from sun exposure and at ambient temperature.

The sterile package should not be stored underneath other packages that could squash it.

Repairs and modifications

TAPES has no parts that can be repaired by the user. TAPES must not be disassembled, modified or repaired by the user. This could hurt patients or damage other medical devices.

Triangulation platform compatibility

TAPES can be used with any triangulation platform with an operating channel measuring 6 to 9 Fr.

It has been validated with Endomina V2.

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