- EN – Instructions for use p.1
- FR – Mode d'emploi p.1
- NL – Gebruikshandleiding p.1
- DE – Bedienungsanleitung p.2
- IT – Istruzioni per l'uso p.2
- ES – Modo de empleo p.3
- PT – Instruções de utilização p.3
- CS – Návod k použití p.4
- SR – Uputstvo za upotrebu p.4
- HU – Használati útmutató p.5

EN – ENGLISH - Instructions for use

Triangulation platform

Do not use Endomina if the package is damaged or unintentionally opened

Terminology and package contents

- A. Arm
- B. Active shoulder
- C. Passive shoulder
- D. Active leg
- E. Deployment leg
- F. Passive leg
- G. Foot
- H. Endoscope fixation wheel
- I. Arm pusher
- J. Brake
- K. Active therapeutic channel
- L. Passive therapeutic channel
- M. Lock connector

Technical information

- Leg length – 860 mm
- Maximum diameter of the part inserted into the patient – 16,5 mm
- Maximum diameter (arm bent, endoscope size not included) – 47 mm
- Shoulder length – 61 mm
- Maximum tool diameter – 1,9 mm
- Maximum duration of use – 4 hours

Intended use

endomina®v2 is a short term, single use medical device for endoscopic gastroplasty using a flexible endoscope and TAPES suturing unit.

Intended users

- The intended users are:
- Nurse(s) trained on the device
- Nurses in training accompanied by clinical field specialist
- Gastro-entérologists trained on the device
- Gastro-entérologists in training accompanied by experienced gastro-entérologist(s).
- All users must understand at least one of the languages in the IFU.

Target patient population

Adults (18 and 65 years)

Indication

Obesity (restrictive surgery on the stomach)

Clinical benefits of endomina used in combination with TAPES

- Excess weight loss
- Improved satiety

Performance characteristics of the device

- endomina®v2 is used for incisionless surgery via the mouth
- endomina®v2 is compatible with TAPES and endoscopes listed in these instructions
- endomina®v2 provides additional working channels
- endomina®v2 offers the ability to perform endoscopic gastroplasty under direct visualization
- endomina offers the possibility to attach/detach from the endoscope within the GI tract without device removal
- surgical intervention with endomina® v2 is completed with short hospitalization.

Contraindications

Women who are pregnant, trying to become pregnant or nursing; children; elderly; pre-existing conditions that may lead to fragile mucosa

For oral insertion:

- Any malformations from mouth to oesophagus (incl. pharynx)
- Any endoscopy contraindication

Warnings

- In order to ensure the safe use of endomina®v2, it is essential to read this paragraph carefully.
- Before using endomina®v2, make sure that the package is properly sealed, make sure that the package has no hole or breach and includes all the components shown in the pictures and that they are in perfect condition. If any component is missing or damaged or if the packaging has been unintentionally opened, do not use endomina®.
- Do not use endomina®v2 after expiration date, it could lead to device malfunction
- If endomina®v2 has been exposed to environmental conditions outside of those specified on the label, do not use endomina®v2.
- endomina®v2 is a device for endoscopic gastroplasty. It can only be used through the mouth
- Always keep at least one additional endomina® in stock before starting a procedure.

- endomina®v2 can only be used by qualified medical staff (a physician or paramedical staff under surveillance by a physician) who has received the necessary training and certification to use it.
- This document describes the use of endomina®v2, not the clinical procedures that can be performed with it.
- endomina®v2 is single-use and cannot be reused.
- endomina®v2 can only be used on patients under general anaesthesia.
- endomina®v2 accepts tools up to 1.9 mm. Do not use larger tools (see "Other endoscopic tools required for the procedure" for more details).
- endomina® can only be used with adult endoscopes (see "Endoscope compatibility" for more details).
- endomina®v2 cannot be resterilized. Resterilizing endomina®v2 can provoke device failure. A resterilisation cycle may not be sufficient to completely eliminate small foreign bodies from hard-to-reach parts (e.g. in long and narrow channels). Moreover, after use, some mechanical components may not fulfil the initial specifications.
- If endomina®v2 is placed outside a clean zone, it cannot be used.
- If one of the components of endomina®v2 becomes detached and falls into a body cavity of the patient, retrieve it with endoscopic forceps.
- Tools can only be inserted into the channels when the distal ends of the shoulders are visible with the endoscope and when the picture is clear, and the arms are straight.
- endomina®v2 may reduce mobility of the endoscope when attached.

Precautions

- Patient should be placed in supine position to avoid entrapment of neighboring organs
- In all cases, an endoscopic examination is necessary prior to using endomina®v2 in order to detect any malformation of the oesophagus.
- If perforated tissue is suspected, it is recommended that an endoscopy be performed to check the state of the tissue and take all necessary measures to treat any lesions prior using endomina®v2.
- Never insert endomina®v2 in the open position.
- Never take endomina®v2 out of the patient in the open position. Similarly, never take endomina®v2 out of the patient without checking visually that endomina®v2 is closed.
- Never open/close endomina®v2 without visual control.
- Always keep the distal ends of the arms in view.
- Never insert or remove endomina®v2 while its arms are bent.
- Always check that endomina®v2 is not too large for the patient's anatomy.
- Never make sudden movements, especially during insertion and removal. This could hurt the patient.
- Always wear appropriate protective clothing, as blood and other tissues carry an infection risk. Suitable protective clothing includes, but not limited to, protective goggles, a mask and waterproof clothes and gloves.
- When using an electric tool (such as a diathermy knife), the other tools must have a nonconductive sheath to avoid electrical shocks for users.
- Do not use force when inserting a tool. If there is resistance, straighten the arm until the tool can be inserted without effort. If this is still not the case, change the tool and/or check that its diameter is small enough for insertion into endomina®v2.
- Do not kink endomina®v2's legs.
- Do not use the endoscope in full retrovision when endomina®v2 is fixed on it.
- Always use CO₂ insufflation.

Potential harms for the patient

- Acute inflammatory tissue reaction
- Bleeding requiring hemostasis management during or after the procedure
- Bleeding (stomach, liver, bowel, gallbladder, pancreas, spleen, omentum)
- Cross contamination
- Delayed intervention
- Infection
- Irritation
- Perforation (esophageal and/or stomach)
- Perforation (pharynx, bowel, bladder, spleen, uterus)
- Perforation (heart, lungs, pancreas, liver, aorta)
- Pneumoperitoneum
- Re-intervention might be required
- Sepsis
- Small mucosal tear
- Strong abdominal pain
- Third party contamination
- User injury

Side effects

Throat irritation

Adverse reactions

The following adverse reactions have been observed during endoscopic operations and may happen during the use of endomina®v2:

- Throat irritation
- Perforation

Procedure

The procedure described here is a check-list for an operation by an experienced user.

1. Check that the endomina®v2 package is not damaged and is perfectly sealed. If any anomaly is detected, do not use endomina®v2; send it back to Endo Tools Therapeutics.
2. For a procedure using oral insertion, perform a prior endoscopic exam to detect any malformation of the oesophagus, which would block endomina®v2.
3. Unpack endomina®v2 and check visually that it is complete and in good condition. Check that no parts protrude.
4. Check that endomina®v2 opens and closes by turning the endoscope fixation wheel (H) and checking visually that the shoulders (B and C) move apart from each other. If they do not move, help the movement by separating them manually. Do not force the movement on the endoscope fixation wheel. Check the action of the arm.
5. Check that endomina®v2's shoulders are closed. If not, close them with the endoscope fixation wheel.
6. Check that the arm is straight. If not, straighten it with the arm pusher.
7. Insert two guide wires into the patient.
8. Insert endomina® (up to 60cm) gently into the patient along the guide wires (in closed position). Use a lubricating gel if needed. Do not use force.
9. Insert the endoscope along endomina®v2's channels.
10. Check visually that there is enough space to open endomina®v2.
11. Open endomina® with the endoscope fixation wheel.
12. Place the endoscope between the endomina®v2 shoulders.
13. Close the endomina®v2 shoulders.
14. Perform the operation according to TAPES IFU.
15. Always straighten the endomina®v2 arm before introducing a new TAPES.

End of procedure

1. Straighten the endomina®v2 arm. If the arm remains bent, see the emergency procedure below.
2. Remove all tools from endomina®v2.
3. Open the endomina®v2 shoulders.
4. Withdraw the endoscope by a few centimetres (2 to 5) to keep the end of endomina®v2 in view.
5. Close the endomina®v2 shoulders. If the shoulders remain open, see the emergency procedure below.
6. Gently remove the endoscope and endomina®v2 from the patient.
7. Dispose of endomina®v2 in an appropriate trash bin.

Emergency procedures

If endomina®v2 gets stuck in the open position inside the patient, use the following procedure:

- If endomina®v2 is stuck on the endoscope:
 1. Turn the endoscope fixation wheel as far as possible (in the opening direction, i.e. unscrewing)
 2. Rotate the endoscope alternately + and – 90° until the endoscope is removed from endomina®v2
 3. If the shoulders can be closed now, continue the normal procedure. Otherwise, continue this procedure.

If endomina®v2 is stuck in the open position, but not on the endoscope:

1. Withdraw endomina®v2 by a few centimeters (2 to 5) to keep the end of endomina®v2 in view.
 2. Turn the endoscope fixation wheel as far as possible (in the closing direction, i.e. screwing in).
 3. If the shoulders are closed enough to pass through the oesophagus (i.e. the space between the shoulders is about 3 mm or less), check that no parts are protruding and gently pull endomina®v2 out, without using force.
 4. If the shoulders are not closed enough to pass through the oesophagus (i.e. the space between the shoulders is still more than about 3 mm), slide the shoulders relative to each other by pushing the passive leg and pulling the handle. This will bend the muscles and bring the shoulders closer to each other. Check that no parts are protruding and gently pull endomina®v2 out, keeping traction between the active and passive legs.
- If endomina®v2 is damaged during the operation and a part stays in the stomach, take it out of the patient.

Storage
endomina®v2 must be stored in its sterile package, in a clean and dry place, away from sun exposure and at ambient temperature. The sterile package should not be stored underneath other packages that could squash it.

Repairs and modifications

endomina®v2 has no parts that can be repaired by the user. endomina®v2 must not be disassembled, modified or repaired by the user. This could hurt patients or damage other medical devices.

Endoscope compatibility

endomina®v2 has been validated with the endoscopes listed at the end of these instructions.

Other endoscopic tools required for the procedure

endomina®v2 can be used with TAPES (Endo Tools Therapeutics).

Devices with metal sheaths are easier to manipulate. It has been validated with the tools listed at the end of these instructions. For other tools, please test prior to procedure. There is no guarantee that instruments selected solely using maximum insertion portion width and working length will be compatible in combination.

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
EN

Any serious incident that has occurred in relation to the medical device should be reported by e-mail to urgent@endotools.be and to the competent authority of the country in which you are established.
Please include your phone number in your e-mail.
Summary of Safety and Clinical Performance is available in the European database on medical devices :
<https://ec.europa.eu/tools/euda/med>
TAPES Basic UDI-DI: 542502928TDS1BZ
IFU is available on www.endotools.be/ifu/

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